Aspirator Suction Unit

Assembly, Installation and Operating Instructions

Yes, you can:

SAVE THESE INSTRUCTIONS

NOTE: Check all parts for shipping damage. In case of damage, DO NOT use. Contact the carrier/Invacare for further instructions.

Safety Summary

To ensure the safe and proper installation of the Aspirator Suction Unit, these instructions **MUST** be followed:

⚠ WARNING

- **DO NOT** install this equipment without first reading and understanding this instruction sheet. If you are unable to understand the warnings, cautions and instructions, contact a healthcare professional, dealer or technical personnel if applicable, before attempting to install this equipment - otherwise, injury or damage may occur.
- Invacare products are specifically designed and manufactured for use in conjunction with Invacare accessories. Accessories designed by other manufacturers have not been tested by Invacare and are not recommended for use with Invacare products.

CAUTION

Federal law restricts this device to sale or rental by or on order of a physician, or any other practitioner licensed by the law of the State in which he/she practices to use or order the use of this device.

⚠ DANGER

- The suction unit is a vacuum suction device designed for the collection of nonflammable fluid materials in medical applications only. Improper use can cause injury or death.
- To reduce the risk of electrocution:
 - ❖ Always unplug this product immediately after using.
 - ❖ **DO NOT** use while bathing.
 - ❖ **DO NOT** place or store product where it can fall or be pulled into a tub or sink.
 - ❖ **DO NOT** place in or drop into water or other liquid.
 - ❖ **DO NOT** reach for a product that has fallen into water. Unplug immediately.

For more information regarding Invacare products, parts, and services, please visit www.invacare.com

⚠ WARNING

- **Proper operation of this product is necessary to avoid physical harm. Improper use can cause personal injury and/or tissue damage. This product should not be used unless the operator has been instructed by a physician or other qualified health care professional.**
 - **DO NOT disassemble this unit. Disassembly or attempted repairs done incorrectly can create an electrical shock hazard. Only qualified technicians should perform servicing of this suction unit. Otherwise, electrical shock or injury may occur.**
 - **To reduce the risk of burns, electrocution, fire or injury to persons:**
 - ❖ **Close supervision is necessary when this product is used by, on, or near children or physically challenged individuals.**
 - ❖ **Use this product only for its intended use as described in this guide.**
 - ❖ **Never operate this product if it has a damaged power cord or plug, if it is not working properly, if it has been dropped or damaged, or if it has been dropped into water. Return the product to a service center for examination and repair.**
 - ❖ **Keep the power cord away from heated surfaces.**
 - ❖ **Never use while drowsy or asleep.**
 - ❖ **Never drop or insert fingers or any other object into any openings.**
 - ❖ **Never block any air openings of the product or place it on a soft surface where the opening may be blocked. Keep all air openings free of lint, dirt, and other foreign objects.**
 - ❖ **DO NOT use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered in a closed environment such as an oxygen tent.**
 - ❖ **This unit is oil-less. DO NOT lubricate.**
 - ❖ **The motor of this compressor has a thermal overload protector. If the motor overheats, the overload protector will shut off the motor. Always disconnect power source before servicing.**
 - ❖ **All electrical products generate heat. To avoid injury, never touch the motor housing during or immediately after operation.**
 - **Some attachments or accessories may not fit the tubing supplied. All attachments or accessories should be checked prior to use to assure proper fit.**
 - **The suction unit is used with a power supply allowing operation on AC voltage (115VAC, 60 Hz). The correct power source must be used to connect to adaptable wall power.**
-

Grounding Instructions

WARNING

DO NOT, under any circumstances, cut or remove the round grounding prong from any plug used with or for Invacare products. Some devices are equipped with three-prong (grounding) plugs for protection against possible shock hazards. Where a two-prong wall receptacle is encountered, it is the personal responsibility and obligation of the customer to contact a qualified electrician and have the two-prong receptacle replaced with a properly grounded three-prong wall receptacle in accordance with the National Electrical Code. If you must use an extension cord, use **ONLY** a three-wire extension cord having the same or higher electrical rating as the device being connected.

Operating Instructions

Set-Up

NOTE: For this procedure, refer to FIGURE 1 on page 4.

1. Place unit on a table top or flat work surface.
2. Check that the power switch is Off (O). Refer to Detail "A".
3. If the unit is cold, let the unit reach room temperature before plugging it in.
4. Connect one end of the connection tubing to the elbow connector on top of the collection bottle, and connect the other end of the connection tubing to the hydrophobic filter.
5. Connect the second length of connection tubing between the hydrophobic filter and the vacuum port barb.
6. Connect the suction tubing to the labeled patient connector (outlet) on the collection bottle.


FIGURE 1 Medical Suction Unit

Using the Suction Unit

CAUTION

Should fluid be aspirated back into the unit, send your suction unit to an authorized service center for a detailed inspection as the vacuum pump may have been damaged.

NOTE: For this procedure, refer to FIGURE 1 and FIGURE 2 on page 5.

NOTE: If the unit sustains a severe drop, validity of the gauge must be checked. Inspect if the float shut-off is still functioning properly when the liquid in the collection bottle reaches the safe full level.

NOTE: Always transport unit with vacuum adjust knob rotated fully clockwise in case unit is dropped.

1. Ensure the suction tubing is connected to the labeled patient connector(outlet) on the collection bottle.
2. Ensure that all connections are secure and without leaks.
3. Connect the unit to power.
4. Turn the power switch On (I).
5. Adjust the vacuum level from 0 to 560 mmHg by turning the vacuum adjust knob:
 - Clockwise to increase vacuum.
 - Counter-clockwise to decrease vacuum.

NOTE: Refer to the vacuum gauge while setting the desired level of vacuum. To accurately read the gauge, block the patient end of the hose or cap off the collection bottle and allow the gauge to reach a stable vacuum reading.

6. Position the suction tubing appropriately with the patient and start therapy.

NOTE: When the fluid reaches the safe full level of the collection bottle (800 cc), the suctioning will stop automatically (float shut-off). Refer to FIGURE 2.

CAUTION

Turn the power switch Off (O) when the automatic float shut-off is activated. Further suctioning could cause damage to the vacuum pump.

The collection jar should be emptied after every use. NEVER let the fluid level exceed the “Safe Full” line. Overflow can damage the aspirator.

7. Turn the power switch Off (O) when the automatic float shut-off is activated.
8. Remove and clean the collection bottle when the motor has completely stopped and the vacuum has dropped. Refer to Cleaning the Collection Bottle and Lid.
9. Unplug the suction unit when not in use.


FIGURE 2 Collection Bottle Components

Cleaning the Collection Bottle and Lid

NOTE: For this procedure, refer to FIGURE 1 on page 4 and FIGURE 2.

NOTE: Clean the collection bottle and lid after each use.

1. Ensure the power switch is Off (O), the motor has completely stopped, and the vacuum has dropped.
2. Disconnect the suction unit from the AC power source.
3. Remove the elbow connector and tubing from the collection bottle lid.
4. Remove the collection bottle lid from the collection bottle.
5. Empty the contents.
6. Clean the collection bottle and lid after each use by one of the following methods:
 - Wash the collection jar and lid in the upper rack of a dishwasher with a water temperature maximum of 150°F (65°C).

-
- Perform the following: (1) Wash with a hot water/dishwashing detergent solution; (2) Rinse with clean, hot tap water; (3) Soak in a one part vinegar to three parts hot water solution for 20 minutes; and (4) Rinse with hot tap water and air dry.
 - Wash with rubbing alcohol and air dry.
 - Wash with a commercial (bacterial-germicidal) disinfectant. Follow disinfectant manufacturer's recommended instructions and dilution rates carefully.

Cleaning the Suction Unit

⚠ WARNING

DO NOT submerge suction unit in water. This will result in severe damage to the vacuum pump.

ALWAYS be sure to remove any foreign material (dirt, lint, etc.) from the suction unit.

1. Disconnect the suction unit from all external power sources.
2. Wipe the suction unit housing with a clean cloth and any commercial (bacterial-germicidal) disinfectant.

Cleaning the Tubing

1. Disconnect the tubing from the connection jar.
2. Clean the tubing after each use by one of the following methods:
 - Wash with rubbing alcohol and air dry.
 - Perform the following (make sure all solutions run through the tubing): (1) Wash in a hot water/dishwashing detergent solution; (2) Rinse with clean, hot tap water; (3) Soak in one part vinegar to three parts hot water solution for 20 minutes; (4) Rinse with hot tap water; and, (5) Air dry.
 - Wash with a commercial (bacterial-germicidal) disinfectant. Follow disinfectant manufacturer's recommended instructions and dilution rates carefully.
3. Keep the outer surface of the tubing clean by wiping with a clean, damp cloth.

Maintenance

⚠ WARNING

Dispose of collection bottle, tubing and hydrophobic filter in a specific waste bin labeled for medical waste. Residues in the collection bottle, tubing and filter may cause respiratory infection.

Inspect tubing and collection bottle for leaks, cracks, etc. before each use.

Replacing the Hydrophobic Filter

NOTE: For this procedure, refer to FIGURE 3.

NOTE: Replace the hydrophobic filter when either of the following conditions occur: (1) the filter is dirty and the aspirator has lost suction, or (2) there is liquid in the hydrophobic filter assembly due to extended use or overflow.

NOTE: The hydrophobic filter should be replaced every two months. However, if overflow occurs, change the filter immediately.

CAUTION

Use only Invacare recommended hydrophobic filters. DO NOT replace with other types of bacteria filters.

1. Remove the elbow connector end of the hydrophobic filter assembly from the top of the collection jar.
2. Separate the opposite end of the hydrophobic filter assembly from the vacuum gauge.
3. Disassemble the hydrophobic filter assembly by disconnecting the two connection tubes from the hydrophobic filter housing.
4. Discard the filter.
5. Re-assemble the hydrophobic filter assembly with a new filter.
6. Re-attach the hydrophobic filter assembly onto the collection bottle lid and vacuum gauge.
7. If necessary, twist the elbow connector to align with the collection jar.


FIGURE 3 Hydrophobic Filter Replacement

Troubleshooting

NOTE: The suction unit contains no user-serviceable parts except the fuse.

PROBLEM	SOLUTION
Suction unit does not turn On.	<ol style="list-style-type: none">1. Check that the plug is firmly fitted to the wall socket.2. Ensure the wall outlet is live. Try plugging in another device such as a lamp.3. Check the power cord for damage.4. Contact your Invacare representative.
Motor runs but no vacuum	<ol style="list-style-type: none">1. Verify tubing connection security.2. Check for leaks or tubing kinks.3. Ensure that float shut-off is not activated.4. Check for bottle leaks and cracks.5. If these solutions do not work, contact your medical provider.
Low vacuum	<ol style="list-style-type: none">1. Use vacuum adjust knob to increase vacuum level.2. Check system for leaks.3. Adjust vacuum adjust knob and release.
Overload protector shuts off the motor frequently	Check number of motor operated appliances on same circuit and adjust for aspirator unit.
Motor overheats and overload protector shuts off motor	<ol style="list-style-type: none">1. Turn power switch Off (O) and allow the motor to cool approximately 5 minutes. Restart compressor.2. Inspect fuses and replace if necessary.3. Call your medical provider if unit still does not start.
Blown fuse	Replace the fuse with the same type of fuse (Miniature Glass Fuse 5x20 mm, fast blow (quick acting), F2A/250V).

Replacement Items/Accessories

PART NO.	DESCRIPTION
IRC1140	Collection jar with 6 ft. suction tubing
1148804	Tube and filter kit

LIMITED WARRANTY

PLEASE NOTE: THE WARRANTY BELOW HAS BEEN DRAFTED TO COMPLY WITH FEDERAL LAW APPLICABLE TO PRODUCTS MANUFACTURED AFTER JULY 4, 1975.

This warranty is extended only to the original purchaser of our products.

This warranty is extended only to the original purchaser who purchases this product when new and unused from Invacare Corporation or a dealer. This warranty is not extended to any other person or entity and is not transferable or assignable to any subsequent purchaser or owner. Coverage under this warranty will end upon any such subsequent sale or other transfer of title to any other person. This warranty gives you specific legal rights and you may also have other legal rights which vary from state to state.

Invacare warrants this product when purchased new and unused to be free from defects in materials and workmanship for a period of one year from date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty. If within such warranty period any such product shall be proven to Invacare's satisfaction to be defective, such product shall be repaired or replaced, at Invacare's option. This warranty does not apply to normal wear and tear or failure to adhere to the product instructions, nor does it include shipping charges incurred in replacement part installation or repair of any such product. Invacare's sole obligation and your exclusive remedy under this warranty shall be limited to such repair and/or replacement.

For warranty service, please contact the Invacare service department during normal business hours at the toll-free telephone number identified on the back cover of this manual. Upon receiving notice of an alleged defect in a product, Invacare will issue a serialized return authorization. It shall be the responsibility of the purchaser to pack the product(s) or part(s) in an appropriate manner to avoid shipping damage and return the product(s) or part(s), at the purchaser's expense, to either Invacare's plant or service center as specified by Invacare in advance. The defective product(s) or component part(s) must be returned for warranty inspection using the serial number as identification within 30 days of the return authorization date. **DO NOT** return any product(s) or part(s) to our factory or service center without prior consent. C.O.D. shipments will be refused - please pre-pay shipping charges.

LIMITATIONS AND EXCLUSIONS: THE FOREGOING WARRANTY SHALL NOT APPLY TO SERIAL NUMBERED PRODUCTS IF THE SERIAL NUMBER HAS BEEN REMOVED OR DEFACED; PRODUCTS SUBJECTED TO NEGLIGENCE, ACCIDENT, IMPROPER OPERATION, MAINTENANCE, OR STORAGE; PRODUCTS MODIFIED WITHOUT INVACARE'S EXPRESS WRITTEN CONSENT (INCLUDING, BUT NOT LIMITED TO, MODIFICATION THROUGH THE USE OF UNAUTHORIZED PARTS OR ATTACHMENTS); PRODUCTS DAMAGED BY REASON OF REPAIRS MADE TO ANY COMPONENT WITHOUT THE SPECIFIC CONSENT OF INVACARE; PRODUCTS DAMAGED BY CIRCUMSTANCES BEYOND INVACARE'S CONTROL; OR PRODUCTS REPAIRED BY ANYONE OTHER THAN INVACARE OR AN AUTHORIZED INVACARE SERVICE CENTER. SUCH EVALUATION SHALL BE SOLELY DETERMINED BY INVACARE.

THE FOREGOING EXPRESS WARRANTY IS EXCLUSIVE AND IN LIEU OF ANY OTHER WARRANTIES WHATSOEVER, WHETHER EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, AND THE SOLE REMEDY FOR VIOLATIONS OF ANY WARRANTY WHATSOEVER, SHALL BE LIMITED TO REPAIR OR REPLACEMENT OF THE DEFECTIVE PRODUCT PURSUANT TO THE TERMS CONTAINED HEREIN. THE APPLICATION OF ANY IMPLIED WARRANTY WHATSOEVER SHALL NOT EXTEND BEYOND THE DURATION OF THE EXPRESS WARRANTY PROVIDED HEREIN. INVACARE SHALL NOT BE LIABLE FOR ANY CONSEQUENTIAL OR INCIDENTAL DAMAGES WHATSOEVER.

SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGE, OR LIMITATION OF HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE EXCLUSION AND LIMITATION MAY NOT BE APPLICABLE.

THIS WARRANTY SHALL BE EXTENDED TO COMPLY WITH STATE/PROVINCIAL LAWS AND REQUIREMENTS.

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All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Invacare Corporation www.invacare.com
© 2007 Invacare Corporation